


GEHPCO Analgesia 350 TENS Unit

**Improve Your
Health Naturally!!!**



GEHPCO

Laguna Hills, California 92653

Email: info@gehpc.com, Tel: (310) 728-6171

GEHPCO

www.gehpc.com
www.gogorders.com
www.alfa1sports.com
www.canteenmwf.com

GEHPCO Analgesia 350 TENS Unit

Analgesia 350 TENS Unit

Model No.: GE-A350

The IntelliScan Pain Management Analog TENS Model PMS349A offers a high-quality, affordable, all-in-one electrotherapy solution. This TENS unit is highly portable, allowing it to be sent home with a patient for their use. It features dual, isolated channels with asymmetrical biphasic output, adjustable pulse and pulse rate, and three modes of operation: conventional, burst or modulation. IntelliScan complete TENS unit will include a 9-volt alkaline battery, two lead wires, electrodes, carrying case and instructions. IntelliScan Pain Management Analog TENS Model PMS349A Features



GEHCPO Analgesia 350 TENS Unit

Analgesia 350 TENS Unit

Model No.: GE-A350



- Features dual, isolated channels with asymmetrical biphasic output, adjustable pulse and pulse rate
- Three-mode operation; conventional, burst or modulation
- Includes 9-volt alkaline battery, two lead wires, electrodes, carrying case and instructions
- Built-in Timer

